

CytoSep™ Cytology Funnels

Single funnels




For SHANDON CYTOSEP cytology centrifuges.

- Used to deposit a thin layer of cells in a clearly defined area of a microscope slide
- Filter paper specially adapted to cytosept, allows the absorption of any excess liquid
- Funnels are supplied with a cap to protect the samples
- For 0.5 ml samples, use single CytoSep™ with white filter paper and cap; 6 mm (28 mm²) deposit - can be used with metal clips 930779
- For samples up to 0.4 ml such as spinal fluid, the simple CytoSep™ with brown filter paper and cap should be used to allow slower absorption of the fluid
- All components are available individually

Cat.No	Description	Pack. unit	Units/ carton	€ Excl. VAT/ carton
930515	Single funnel with white filter and cap	50	500	NC -
930663	Single funnel with white filter and cap, individually packed	100	500	NC -
930774	Single funnel with brown filter and cap	50	500	NC -
930775	Single funnel only	-	500	NC -
930776	Cap only	50	500	NC -
930596	White filter paper	-	200	NC -
930784	Brown filter paper	-	200	NC -

Double funnels



For 2 samples simultaneously on one slide.

- Ideal for immunohistochemistry
- For samples up to 0.5 ml; 6 mm deposit (28 mm²)
- The filter is pre-attached
- Delivered with cap
- Can be used with metal clips 930779
- All components are available individually

Cat.No	Description	Pack. unit	Units/ carton	€ Excl. VAT/ carton
930664	Double funnel with white filter and cap	50	500	NC -
930665	Double funnel with white filter and cap, individually packed	100	500	NC -
930777	Double funnel only	-	500	NC -
930776	Cap only	50	500	NC -
930778	White filter paper	-	200	NC -

Metal clip for CytoSep™

Keeps the slide close to the CytoSep™ during centrifugation in the CytoSpin™

- Autoclavable

Cat.No	Description	Units/ carton	€ Excl. VAT/ carton
930779	Metal clip	6	NC -



CytoSep™ Single or double funnels made entirely of plastic



- Single use
- Reduce the risk of exposure to pathogens and cross-contamination
- Do not require the use of metal clips
- For 0.5 ml samples, use the single Cytosept with white filter paper and cap; 6 mm deposit (28 mm²)
- Double funnels allow the deposition of 2 samples simultaneously on one slide
- Supplied with filter paper
- All components are available individually

Cat. No.	Description	Pack. unit	Units/ carton	€ Excl. VAT/ carton
930543	Single funnel with white filter and cap	40	480	NC -
930780	Single funnel with brown filter and cap	40	480	NC -
930544	Double funnel with white filter and cap	40	480	NC -
930781	Single funnel only	50	500	NC -
930782	Double funnel only	50	500	NC -
930776	Cap only	50	500	NC -
930783	Base only	50	500	NC -
930596	White filter paper for single funnel	-	200	NC -
930784	Brown filter paper for single funnel	-	200	NC -
930778	White filter paper for double funnel	-	200	NC -



ANALYSIS AND MEASUREMENTS