


CytoSep™ Cytology Funnels (continued)

CytoSep™ Mega funnel



For the preparation of high volume samples.


- Does not require the use of metal clips
- Large deposition area : 22 x 14,75 mm (325 mm²) ; 12 times the volume of the standard CytoSep™ (6 ml)
- Strategically placed baffles prevent cell sedimentation, resulting in uniform cell deposition and excellent slide preparation quality
- Reduce the risk of exposure to pathogens and cross-contamination

| Cat. No. | Description | Pack. unit | Units/ carton | € Excl. VAT/ carton |
|----------|------------------------------|------------|---------------|---------------------|
| 930785 | Cytosep Mega funnel with cap | 40 | 480 | NC - |
| 930786 | White filter paper | - | 200 | NC - |

CytoSep™ for Unitech (Wescor) Cytopro Centrifuge

➤ Single funnels

- The individual elements fit together quickly and are easily aligned
- Compression ring for better flow control and results
- Samples are loaded directly through the cap
- The caps provide additional safety for the operator
- Large, centred cell deposition area for easier and more sensitive screening



| Cat. No. | Description | Pack. unit | Units/ carton | € Excl. VAT/ carton |
|----------|--|------------|---------------|---------------------|
| 930791 | Single chamber with white filter paper | 24 | 48 | NC - |
| 930792 | White filter paper | - | 100 | NC - |

➤ Double funnels

- Double chamber with close and easy to locate deposit areas
- 7 mm spotlight provides a 37% larger area of cell collection



| Cat. No. | Description | Pack. unit | Units/ carton | € Excl. VAT/ carton |
|----------|--|------------|---------------|---------------------|
| 930793 | Double chamber with white filter paper | 24 | 48 | NC - |
| 930794 | White filter paper | - | 100 | NC - |

CytoSep™ funnels for Sakura Cyto-Tek centrifuge



| Cat. No. | Description | Pack. unit | Units/ carton | € Excl. VAT/ carton |
|----------|--------------------|------------|---------------|---------------------|
| 930787 | Chamber 1 ml | 50 | 200 | NC - |
| 930788 | Support | 50 | 200 | NC - |
| 930789 | Cap | 50 | 200 | NC - |
| 930790 | White filter paper | - | 200 | NC - |

CytoSep™ funnels for Cytofuge StatSpin centrifuge



- For gentle concentration of cells in a 7 mm diameter area on a slide
- Supernatant absorbed by the filter
- For samples from 50 to 400 µl
- High quality single cell coating
- Filter sold separately

| Cat. No. | Description | Units/ carton | € Excl. VAT/ carton |
|----------|--------------|---------------|---------------------|
| 930795 | Concentrator | 48 | NC - |
| 930796 | Filter | 200 | NC - |

CytoSep™ funnels for Hettich Cyto-System



- Multiple funnel options: 2, 3 or 4 per slide
- 4 sizes available : 1, 2, 4 or 8 ml
- The 1 ml chamber is intended for the small volume with low cell concentration
- 2 and 4 ml chambers: suitable for cell-rich samples
- 8 ml chamber: recommended for very large samples (urine)

| Cat. No. | Description | Units/ carton | € Excl. VAT/ carton |
|----------|-----------------------------|---------------|---------------------|
| 930797 | Funnel chamber 1 ml | 50 | NC - |
| 930798 | Funnel chamber 2 ml | 50 | NC - |
| 930545 | Funnel chamber 4 ml | 50 | NC - |
| 930799 | Funnel chamber 8 ml | 50 | NC - |
| 930800 | Chamber filter 1.2 and 4 ml | 200 | NC - |
| 930801 | Chamber filter 8 ml | 200 | NC - |