



# TRANSFER PIPETTES

Select your transfer pipette according to the volume of the bulb required.

	FINE TIP		STRAIGHT ROD				SPECIAL		
Bulb Volume (ml)	4.8	3.5	0.9	3.3	8	10	0.2	0.2	0.2
<b>Length (mm)</b>	<b>147</b>	<b>157</b>	<b>64</b>	<b>155</b>	<b>155</b>	<b>171</b>	<b>131</b>	<b>131</b>	<b>124</b>
Drop per ml	45	50	25	25	19	18	-	-	-
Bulb Volume + stem (ml)	8.7	5.8	1.2	4.6	15	-	1	1	0.8
Graduation (ml)	-	-	-	-	-	-	-	-	-



 Drop volume (µl)

<b>Non-sterile</b>	<b>043234</b>	<b>043235</b>	<b>043241</b>	<b>043251</b>	<b>043252</b>	<b>030126</b>	<b>043292</b>	<b>043293</b>	<b>043295</b>
Box of	400	500	500	500	250	200	500	500	500
€ Excl. VAT/ box	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -
<b>Sterile/1</b>	<b>043234S</b>	<b>043235S</b>	<b>043241S</b>	<b>043251S</b>	<b>043252BS</b>	-	-	-	-
Box of	400	400	500	500	250	-	-	-	-
€ Excl. VAT/ box	NC -	NC -	NC -	NC -	NC -	-	-	-	-
<b>Sterile/10</b>	<b>043234M</b>	<b>043235M</b>	<b>043241M</b>	<b>043251M</b>	<b>043252M</b>	-	-	-	-
Box of	400	500	500	500	250	-	-	-	-
€ Excl. VAT/ box	NC -	NC -	NC -	NC -	NC -	-	-	-	-

ANALYSIS AND MEASUREMENTS