

# INDICATORS FOR STERILISATION


## Autoclave indicator tape



- Self-adhesive
- Roll 19 mm x 50 m
- For steam sterilisation between 121 °C and 134 °C
- The tape changes colour when sterilisation is carried out

Cat. No.	Description	€
490009	Autoclave indicator tape	NC -


## Neutral autoclave adhesive tape



Adhesive tape for closing containers or packaging products for autoclaving. Does not have a sterilisation indicator.

Cat. No.	Description	€
150262	18 mm x 50 m roll	NC -
150272	24 mm x 50 m roll	NC -
150263	48 mm x 50 m roll	NC -
086000	Tape dispenser	NC -

## Sterilisation indicator tapes



The adhesive indicator tapes to be placed on the products to be sterilised allow the passage to sterilisation to be certified. Coloured bands appear at the end of the sterilisation process. 3 models of indicator tapes are available: for autoclave (wet heat), for poupinel (dry heat) and for gas (ethylene oxide).

Dispensers 086000 and 140265 are suitable for use with the tape (see page 1235).

Cat. No.	Description	Dimensions	Units/box	€/box
150260	Autoclave Indicator tape	50 m x 19 mm	1	NC -
150261	Poupinel Indicator tape	55 m x 19 mm	1	NC -
095038	Ethylene oxide Indicator tape	50 m x 19 mm	1	NC -
150267	Autoclave Indicator tape 3M	55 m x 18 mm	1	NC -
150283	Autoclave Indicator tape 3M	55 m x 24 mm	1	NC -

## Steripoint® sterilisation indicator



- High accuracy in the shift time
- Constant reproducibility
- Quality commitment due to the numerous controls of shift times during production
- The shift of Steripoint® is inhibited by the presence of a steam / air or gas mixture, or when one of the parameters of the cycle is insufficient
- The shift from purple to green ink indicates compliance of time, temperature and steam conditions during sterilisation cycles

Cat. No.	Colours	Sterilisation point	Units/box	€/box
950171	Blue	134 °C 3,5 min. - 121 °C 15 min.	250	NC -
950172	Pink	134 °C 5,3 min. - 121 °C 15 min.	250	NC -
950173	Orange	134 °C 7 min. - 121 °C 20 min.	250	NC -
950174	Green	134 °C 18 min.	250	NC -

## Sterilisation and washing indicators



- Non-sequential integrator
- Declination according to autoclave cycles :
  - TST SKC short cycle : 121 °C/15 min
  - TST 3702 standard cycle : 125 °C/20 min

Cat. No.	Description	Units/box	€/box
972446	Sterilisation and washing indicator TST SKC short cycle	200	NC -
972440	Sterilisation and washing indicator 3702 SKC short cycle	200	NC -