

## Steranios 2% and Opaster'Anios - Total cold disinfection of instruments



- ▶ **Instructions for use**
  - Ready to use: do not dilute
  - Clean and rinse all instruments very thoroughly before soaking
  - Soaking time according to the desired efficiency spectrum:
    - bactericidal, fungicidal, virucidal : 20 minutes
    - sporicidia : 1 hour

Cold disinfectant for instruments.

### ▶ Steranios 2%

- ▶ **Simple: no neutral pH activator**
- ▶ **Effective: bactericidal, fungicidal and virucidal in 20 minutes, sporicidal in 1 hour**

### ▶ Opaster'Anios

- ▶ **High level disinfectant for thermosensitive medical devices**
- ▶ **Bactericidal, virucidal, fungicidal, and mycobactericidal activity in 5 minutes**
- ▶ **Anti-corrosive effect on metals sterile or filtered 0.2 µm**
- ▶ **Wide compatibility with all types of materials**
- ▶ **Active on Clostridium Difficile 0127**

### Steranios 2%

Cat. No.	Description	€
059604	Steranios 2%, 5 litre canister	NC -
059605	Steranios 2%, 4 x 5 litre canister	NC -

### Opaster'Anios

Cat. No.	Description	€
059917A	Opaster'Anios, 5 litre canister	NC -
059917	Opaster'Anios, 4 x 5 litre canisters	NC -
059987	Control strips - tube of 50	NC -

## ANIOS soaking bins



Soaking containers for collection, pre-disinfection of instrumentation.

- ▶ **Two models: 2 and 5 litres**
- ▶ **For use with SterAnios and Aniosyme DD1**
- ▶ **2 litre version with or without slot on the lid**

Cat. No.	Model	Int. dim. (mm)	Ext. dim. (mm)	€
059815	2 l.	230 x 170 x 90	295 x 200 x 130	NC -
059817	2 l. with slot	230 x 170 x 90	295 x 200 x 130	NC -
059816	5 l.	365 x 255 x 125	430 x 290 x 125	NC -



Find all our soaking containers on [www.dulis.nl](http://www.dulis.nl)

## Autoclavable soaking container



- ▶ **For pre-disinfection and disinfection**
- ▶ **Polypropylene tray autoclavable at 134°C for 18 minutes**
- ▶ **Equipped with a fish bowl and a lid**

Cat. No.	Description	Dimensions L x W x H (mm)	Volume (l)	€
972534	Tray	257 x 180 x 128	3	NC -
972536	Bin + tray	890 x 200 x 145	12	NC -

## Non-autoclavable soaking container



- ▶ **For pre-disinfection and disinfection**
- ▶ **Equipped with a fish bowl and hollow handles**
- ▶ **Transparent anti-splash lid for easy insertion of objects without opening the lid**

Cat. No.	Description	Dimensions L x W x H (mm)	Volume (l)	€
972538	Tray	550 x 200 x 120	5	NC -
972539	Container without tap	390 x 290 x 170	10	NC -
972540	Container with tap	615 x 400 x 220	30	NC -