

CORNING

Corning® Spinner



Spinner with deflectors and stirring blades for optimal aeration and mixing of the culture.

- Ready to use, supplied with stirring paddle and magnet bar
- Sterile polystyrene (gamma radiation)
- Individual packaging
- Easy access for serological pipettes

Wide choice of volumes and materials, consult us



Cat. No.	Total volume (ml)	Ø x h (mm)	Material	Cap Type	Units/bag	Units/carton	€ Excl. VAT/carton
003152	125	63,5 x 145	PS	Standard	1	12	NC -
003153	500	87,3 x 203,2	PS	Standard	1	12	NC -
003561	1000	137 x 245	PS	Standard	1	6	NC -
003563	3000	188 x 259	PS	Standard	1	4	NC -
003578	500	63,5 x 145	PS	Filter	1	12	NC -
003580	1000	137 x 245	PS	Filter	1	6	NC -
003581	3000	188 x 259	PS	Filter	1	4	NC -



Duran® GLS80 reactors with magnetic stirring rod



091909



091912

GL14 and GL18 fittings for the introduction and removal of materials during mixing, for the connection of external diameter 1.6 - 12 mm or sterile pressure compensation.

- Much better mix than standard magnet bars
- The stirring rod is interchangeable and can be used at a speed of 500 rpm
- Maximum operating temperature of 140 °C
- Autoclavable
- Easy access for serological pipettes

Applications

- Mixing of liquids
- Dissolution of solids
- Simple fermentation process

Accessories

Cat.No.	Désignation	€ Excl. VAT
091912	Helix magnet for GLS80 cap	NC -
091913	Spare anchor magnet for GLS80 cap	NC -
091914	Spare GLS80 cap for stirring rod	NC -
091915	Spare stainless steel shaft with PEEK thread for stirring rod for GLS80 cap	NC -

Cat.No.	Désignation	€ Excl. VAT
091911	Plug GLS80 with magnetic anchor stirring rod	NC -
096328	Duran® 250 ml bottle with PP GLS80 blue/gray cap, 1 PP GL14 blue cap, 2 GL14 PBT red caps, 2 PBT GL18 red caps	NC -
091909	Duran® 1 L bottle with PP GLS80 blue/gray cap, 1 PP GL14 blue cap, 2 GL14 PBT red caps, 2 PBT GL18 red caps	NC -
091910	Duran® 2 L bottle with PP GLS80 blue/gray cap, 1 PP GL14 blue cap, 2 GL14 PBT red caps, 2 PBT GL18 red caps	NC -



Duran® GL45 reactors with magnetic stirring rod



GL14 and GL18 fittings for the introduction and removal of materials during mixing, for the 1.6 - 6 mm outer diameter connection or sterile pressure compensation

- Much better mix than standard magnet bars
- Can be used at a speed of 500 rpm
- Maximum operating temperature of 140 °C
- Autoclavable

Applications

- Mixing of liquids
- Dissolution of solids

Cat. No.	Désignation	€ Excl. VAT
091982	Duran bottle® 500 ml with PP GL45 blue/gray cap, 1 PP GL14 blue cap, 2 GL14 PBT red caps	NC -
091983	Duran Bottle® 1000 ml with PP GL45 blue/gray cap, 1 PP GL14 blue cap, 2 GL14 PBT red caps	NC -

Accessories

Cat.No.	Désignation	€ Excl. VAT
091985	Retractable eyelet magnet for GL45 cap	NC -
091984	Spare GL45 blue/grey cap with PP GL14 blue cap	NC -