




Molecular biology grade agarose powder



- > Ultra pure agarose with low background noise
- > Easy to dissolve
- > Can be used in low concentrations (0.75 - 2%)
- > Ideal for the separation of fragments from 250 bp to 25 kb
- > Gel strength : < 1000 g/cm (1 %), > 2000 g/cm (1.5 %)
- > Clarity : < 4 NTU



Cat. No.	Description	€ Excl. VAT
523277	Agarose 100 g	NC -
523278	Agarose 500 g	NC -



Enter the Cat. No. on www.dulis.nl to find the Agarose selection guide



Basic agarose

- > Used for the separation of DNA and RNA fragments at concentrations ranging from 0.8 to 2%

Cat. No.	Description	€ Excl. VAT
348376	Agarose 500 g	NC -

Agarose



Cat. No.	Type	Applications	Flask	€ Excl. VAT
4905008	High.melting temp./medium.resol..	Standard on DNA	25 g	NC -
4905009	High. melting temp./medium. resol..	Standard on DNA	500 g	NC -

