

Multiflex and Miniflex capillary tips for gel



Fine capillary tips for high precision gel deposition.

- Round-ended tips (ø ext. : 0.57 mm) are particularly suitable for SDS PAGE gels.
- The spatula tips with flat end (ø ext. : 0.17 and 0.37 mm) are indicated for the sequencing gel. The great length tips 200 µl tips (83 mm) allow for use in deep tubes.

► Loaded on an autoclavable rack with a lid for 200 tips



For all brands of 200 µl pipettes

Volume (µl)	Tip	Out. ø (mm)	Length (mm)	On 200-tip rack Non-sterile			On 200-tip rack Sterile		
				Cat. No.	Units/box	€ Excl. VAT/box	Cat. No.	Units/box	€ Excl. VAT/box
0.5 - 200	Rounded	0.57	83.5	021570	2 X 200	NC -	013810	4 X 200	NC -
0.5 - 200	Spatulated flat	0.37	83.5	027600	2 X 200	NC -	013770	4 X 200	NC -
0.5 - 200	Spatulated flat	0.17	83.5	017310	4 X 200	NC -	017320	4 X 200	NC -



For 10 µl micropipettes Gilson, Sartorius Biohitfamily, Labsystems

Volume (µl)	Tip	Out. ø (mm)	Length (mm)	On 200-tip rack Non-sterile			On 200-tip rack Sterile		
				Cat. No.	Units/box	€ Excl. VAT/box	Cat. No.	Units/box	€ Excl. VAT/box
0.1 - 10	Rounded	0.57	45	015100	4 X 200	NC -	015110	4 X 200	NC -
0.1 - 10	Spatulated flat	0.37	45	015120	4 X 200	NC -	015130	4 X 200	NC -
0.1 - 10	Spatulated flat	0.17	45	017350	4 X 200	NC -	017360	4 X 200	NC -



For 10 µl Eppendorf micropipettes

Volume (µl)	Tip	Out. ø (mm)	Length (mm)	On 200-tip rack Non-sterile			On 200-tip rack Sterile		
				Cat. No.	Units/box	€ Excl. VAT/box	Cat. No.	Units/box	€ Excl. VAT/box
0.1 - 10	Spatulated flat	0.37	83.5	-	-	-	013900	4 X 200	NC -
0.1 - 10	Spatulated flat	0.17	83.5	017330	4 X 200	NC -	017340	4 X 200	NC -



Multiflex capillary tip 0.5 - 200 µl bulk

Cat. No.	Description	Units/box	€ Excl. VAT/box
028480A	Multiflex Tip 0.5 - 200 µl	1000	NC -

Extra long 10 µl universal tip

Cat. No.	Description	Units/box	€ Excl. VAT/box
077525	Tip, in bulk	1000	NC -

Long reach tip (44 mm), for Gilson, Eppendorf, Sartorius Biohitfamily, Labsystems. Ideal for micro-volumes in universal molecular biology.

