

EpMotion® 96 and 96XL

- Semi-automated pipettor for all types of 96 and 384 well plates
- Automatic recognition of EpTIPS Motion tips available with or without filter
- Minimum space requirement
- Downloadable programming and control software
- 10 predefined and configurable methods

- **Standard model** : 0.5 µl to 300 µl available in 1 or 2 positions
- **Model XL** : 5 µl to 1000 µl available in 1 or 2 positions

	96		96XL	
Number of positions	1	2	1	2
Cat. No.	999506	933947	933948	933949
€ Excl. VAT	NC -	NC -	NC -	NC -

EpTIPS Motion tips

Cat. No.	Volume (µl)	Sterile	Filter	PCR clean	Units/box	€ Excl. VAT/box
999513	1 - 50	No	Yes	Yes	24 x 96	NC -
999516	20 - 300	No	Yes	Yes	24 x 96	NC -
999514	20 - 300	No	No	No	24 x 96	NC -
933223	40 - 1000	No	Yes	Yes	24 x 96	NC -
933222	40 - 1000	No	No	No	24 x 96	NC -

eppendorf



Accessories

Cat. No.	Description	€ Excl. VAT
933951	195 ml PCR Clean tank with low dead volume	NC -
933952	2-position slider upgrade kit for EpMotion 96	NC -

ZipTip pipette tips

For the concentration and purification of biological samples

The ZipTip unit is a 10 µl pipette tip with a bed of chromatographic media of approximately 0.6 µl attached to the end. It is ideal for concentrating, desalting and fractionating protein, peptide and oligonucleotide samples prior to mass spectrometry, HPLC, capillary electrophoresis and other analytical techniques.

Several resins are available depending on your application:

- C18/µ-C18 :
 - desalting and concentration of peptides, proteins or oligonucleotides
 - fractionation of complex mixtures of proteins or peptides
- C4: desalting and protein concentration.
- SCX - cation exchange :
 - purification of peptides
 - removal of detergents

Box with	8	96	960
C18	051168	051177	051182
µ-C18	051169	051178	051275
C4	051170	051465	051183
€ Excl. VAT/box	NC -	NC -	NC -
SCX	051468	051368	-
€ Excl. VAT/box	NC -	NC -	-

