


AlumaSeal® 96 et AlumaSeal® 384 ClearLine® Films

ClearLine®

- Applications : PCR, long-term storage
- 38.1 µm aluminium
- Operating temperature : -80 °C to +120 °C
- With a tab
- Drillable
- Opaque
- Reinforced sealing for 96-well plates



Cat. No.	Description	Dimensions (mm)	Units/ carton	€ Excl. VAT/ carton
760231	AlumaSeal® 96 ClearLine® Film for non-skirted 96 plate	117.5 x 77.8	100	NC -
760232	AlumaSeal® 384 ClearLine® Film for 384-well or deep-well plate	123.8 x 82.6	100	NC -




MOLECULAR BIOLOGY

Alumaseal® CS™ Adhesive film for freezing

ClearLine®


- For cold storage of microplates
- Dimensions : 81.5 x 132.6 mm
- Length without tabs : 123.04 mm
- Resistant to -80°C

Cat. No.	Description	Units/ bag	Units/ carton	€ Excl. VAT/ carton
760243	Film Alumaseal® CS™	100	100	NC -
760244	Film Alumaseal® CS™ sterile	25	50	NC -



FoilSeal™ ClearLine® Film

ClearLine®



- Applications : PCR
- 50.8 µm aluminium
- Operating temperature : -20 °C to 120 °C
- Universal
- Dimensions without tab: 118.1 x 78.3 mm
- With a detachable tab
- Drillable
- Opaque

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
765282	Film FoilSeal™ ClearLine®	100	NC -
765281	Film FoilSeal™ ClearLine® sterile	50	NC -

EZcap™ Film 1 x 8 PCR FoilSeal™ ClearLine®

ClearLine®



- Applications : PCR
- 50.8 µm aluminium
- Operating temperature : -20 °C to 120 °C
- For 8-tube strips and plates
- Dimensions without tab: 81.2 x 10.1 mm
- With two detachable tabs
- Drillable
- Opaque

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
765268	EZcap™ Film 1 x 8 PCR FoilSeal™ ClearLine®	300	NC -



eXTReme™ FoilSeal™ ClearLine® Film for PCR

ClearLine®

- Universal aluminium waterproof film
- Thickness : 50 µm
- Dimensions : 78.3 x 135.9 mm
- Length without tabs : 117.1 mm
- Excellent chemical resistance (including DMSO)
- Temperature : -20 °C to +120 °C
- For PCR, long term storage, light sensitive experiments and robotic applications
- Drillable with pipette tips and robotic tips
- Removes without leaving a residue
- Single film with one detachable tab, roll film with two detachable placing tabs (width 10 mm)



Cat. No.	Description	Dimensions (mm)	Units/ carton	€ Excl. VAT/ carton
765266	eXTReme™ FoilSeal™ Film for PCR	78,3 X 129	100	NC -
765270	eXTReme™ FoilSeal™ Film roll for PCR	78,3 X 139	500	NC -