

Bottles for specific applications

Borosilicate glass carboys with tap


Graduated borosilicate glass carboys with tubing at the base of the GL32 thread.

Screw cap neck DIN GL45. Supplied complete with screw connection cap, silicone seal, flow tap and PTFE needle valve.

Cat. No.	Volume (l)	€
046420	1	NC -
046421	2	NC -
046422	5	NC -
046423	10	NC -

Replacement taps

Cat. No.	Description	€
046424	Spare tap for 1 and 2 litre carboys	NC -
046425	Spare tap for 5 and 10 litre carboys	NC -



Tilt Duran® Cell Culture Flasks

- ▶ Borosilicate glass 3.3
- ▶ GL56 white polypropylene cap, Ø 56 mm, autoclavable
- ▶ 2 positions: upright (filter sterilization and storage) and tilted at 45° (pipetting)
- ▶ Graduations visible on both sides and in both positions
- ▶ DURAN® Tilt Light shield, white silicone sleeve
- ▶ Identification rings : GL56, 4 colours, autoclavable
- ▶ Vacuum filtration unit for Tilt with adapter GL45 (ext)/GL56 (int), 507 mbar recommended, 677 mbar max.



Cat. No.	Description	Height (mm)	Volume (ml)	Units/box	€/box
096014	Tilt bottle with white GL56 cap (non-sterile)	151	500	4	NC -

Accessories

Cat. No.	Description	Material	Colours	Units/box	€/box
096015	Protective shell, 4 identification rings included	Silicone	Orange, yellow, blue, pink	4	NC -
096016	Identification ring GL56	Silicone	Orange	20	NC -
096017	Identification ring GL56	Silicone	Yellow	20	NC -
096018	Identification ring GL56	Silicone	Blue	20	NC -
096019	Identification ring GL56	Silicone	Pink	20	NC -
096020	Adapter GL45 for 45 mm filter unit (autoclavable)	PTFE	White	1	NC -
096021	Self-adhesive label for GL56 cap, Ø 40 mm	Polyester	White	60	NC -
096022	GL56 Cap	Polypropylene	White	10	NC -

Find photos of the accessories on www.dulis.nl



TILT Bottle Top Filter Unit

- ▶ Capacity : 500 ml
- ▶ For GL45 neck bottles
- ▶ Gamma sterilised
- ▶ PES membrane, 90 mm
- ▶ Different porosities available

Cat. No.	Porosity (µm)	Units/box	€/box
096023	0.1	12	NC -
096024	0.2	12	NC -
096025	0.45	12	NC -

Pyrex® GL 56 bottle, dual position

- ▶ Borosilicate glass
- ▶ Can be placed in a vertical position or in 45° angle position for ease of use
- ▶ Bottle is graduated to 500ml in the vertical position and 400ml in 45° angle position
- ▶ Trace Code
- ▶ Can be sterilised chemically or by dry heat

Cat. No.	Capacity (ml)	Bottle height (mm)	Thread size (mm)	Units/box	€/box
785332	500	151	GL 56	4	NC -

