

MORTARS AND PESTLES


Agate Mortars, Pestles and Beads



Mortar and pestle

Cat. No.	Inner ø (mm)	Outer ø (mm)	Depth (mm)	Volume (ml)	€
140870	35	40	10	5	NC -
140871	82	100	25	70	NC -
140872	100	120	30	140	NC -
140873	110	130	33	150	NC -
140874	125	150	38	200	NC -
140875	128	160	39	275	NC -

Agate Mortars and Pestles



Cat. No.	Inner ø (mm)	Depth (mm)	Tot. Height (mm)	Capacity (ml)	Pestle length (mm)	€
068543B	50	14	26.5	15	59	NC -
068544	61	20	32	29	72	NC -
068545B	63	23	38	30	83	NC -

Porcelain Mortars and Pestles



Mortars only			Smooth interior		Rough interior	
Capacity (ml)	Inner ø (mm)	Int. height (mm)	Cat. No.	€	Cat. No.	€
30	56	36	140619	NC -	140631	NC -
70	63	41	140620	NC -	140632	NC -
75	76	45	140621	NC -	140633	NC -
110	83	50	140622	NC -	140634	NC -
160	90	56	140623	NC -	140635	NC -
220	105	64	140624	NC -	140636	NC -
400	125	71	140625	NC -	140637	NC -
500	139	82	140626	NC -	140638	NC -
650	150	90	140627	NC -	140639	NC -
1000	180	104	140628	NC -	140640	NC -

Single pestles for porcelain mortars

Cat. No.	ø (mm)	Length (mm)	For mortars ø (mm)	€
140643B	23	115	70, 81 and 91	NC -
140644B	26	135	100 y 113	NC -
140645B	37	150	125 y 144	NC -
140646B	44	170	163	NC -
140647	58	210	183	NC -
140648	73	255	215	NC -
140649	86	305	257	NC -

Cryogenic mortar



Cryogenic mortar version

- Supplied without pestle
- Dimensions (ø x H): 16.5 x 11.4 cm
- Handle for removing stainless steel mortar
- Blue polyethylene bowl containing liquid nitrogen
- Suction cups for mortar stability

Mini Cryogenic mortar version

- Allows the grinding of small frozen tissue samples contained in a 1.5 ml microtube
- Supplied with a stainless steel pestle with a polyethylene handle and a handle for removing the mortar

Cat.No.	Description	€
088763	Cryogenic mortar	NC -
088767	Cryogenic Mini mortar	NC -