## TissueMaster-125 Omni Homogeniser



- ➤ Power : 125 W
- ➤ Economical for samples from 0.25 to 30 ml
- ➤ Variable speed motor from 5,000 to 35,000 rpm
- ➤ Supplied with a 7 mm x 85 mm C189 dispersion rod
- Most samples are dispersed in less than a minute
- ➤ Easy to clean and dismantle
- ➤ Portable and lightweight
- ➤To be handled by hand or mounted on a stand
- >CSA, UL and CE C196 standards
- ➤ Noise below 72 dB
- >1 year warranty for the motor

Cat. No.	Description	€ Excl. VAT
078053	TissueMaster-125 Homogeniser withrod 7 x 85 mm	NC -

## **TH Omni Homogeniser**



- ➤ To be held in the hand or on a stand
- > For samples from 0.05 to 100 ml
- ➤ Lightweight, ergonomic
- >Supplied with probe/dispersion rod assembly tool and manual
- ▶125 W motor
- ➤ Mounting of EZ connect stainless steel probes with 1/4 revolution bayonet
- ➤ Noise limited to 72 dB
- ➤ Variable speed from 5,000 to 35,000 rpm
- **▶**CE and CSA certified

TH pack with TH 220 V homogeniser motor , 3 EZ-Connect rods for fine grinding (5 x 75 mm flat bottom (Cat No. 078007), 7 x 115 mm saw teeth (Cat. No. 078148) and 10 x 115 mm saw teeth (Cat. No.. 078150), nut and stand.

-			
Cat. No.	Description	€ Excl. VAT	
078000	Full TH pack	NC -	
078001	Omni Tissue Homogeniser only	NC -	
078002	Stand	NC -	
078003	Rod with tightening nut	NC -	

## **THQ Omni Homogeniser**



- ➤ For samples from 200 µl to 120 ml
- ➤To be held by hand or on a stand
- ➤ Variable speed from 8,000 to 32,000 rpm
- Max. noise level 68 dB
- ▶40 W motor
- > Supplied with assembly tool, stand clamp, manual and starter kit including 12 plastic dispersion rods, 6 for soft tissue and 6 for hard or frozen tissue
- ➤ Dispersion rod :
  - Ø : 7 x 110 mm
  - Recommended working volume : o.25 A 30 ml

## Omni GLH-850 Homogeniser



- For samples from 1.5 ml to 5 l
- >850 W motor, 2 year warranty
- ➤ Variable speed : 5000 to 25000 rpm
- ➤ Noise level : < 72 dB
- ➤Timer up to 9 min 59 sec
- ➤On / off button
- ➤Base included
- ➤CE certified

Cat. No.	Description	€ Excl. VAT
078211	THQ homogeniser	NC -
078002	Stand	NC -

CatNo.	Description	€ Excl. VAT
	GLH-850 homogeniser complete with 3 stainless steel rods (10, 20 and 30 mm)	NC -
078147	GLH-850 Single motor homogeniser	NC -