

Mini Cooler True North

- Cooling boxes with insulating gel
- Temperature maintenance or rapid cooling of samples
- Polycarbonate with non-slip feet
- Temperature maintenance graph delivered with each version
- Lid with marking grid



Cat. No.	Temperature	Number of tubes	Tube type	Lid	L x W x H (mm)	Temperature maintenance	Colours	€ Excl. VAT
320538	0 °C	3 x 4	0,5* and 1,5 - 2 ml microtubes or cryotubes	Transparent	145 x 105 x 100	1 °C for 4 hours.	Green	NC -
320539		4 x 8		Transparent	245 x 105 x 100	1 °C for 5 hours.		NC -
320540		4 x 8		Insulation	245 x 105 x 100	1 °C for 7,5 hours.		NC -
320541	- 20 °C	3 x 4		Transparent	145 x 105 x 100	- 15 °C for 2,5 hours.	Blue	NC -
320542		4 x 8		Transparent	245 x 105 x 100	- 15 °C for 3,5 hours.		NC -
320543		4 x 8		Insulation	245 x 105 x 100	- 15 °C for 5 hours.		NC -
320544		8 x 12	0.2 ml microtubes or 96-well plate	Transparent	225 x 165 x 100	- 15 °C for 3,5 hours.		NC -

* with adapters supplied

Nalgene Labtop coolers

Thermo
SCIENTIFIC | NALGENE

- Maintain the temperature of samples on the bench at 0°C or -20°C
- In polycarbonate
- Accept microtubes or tubes depending on the model
- The case is filled with a non-toxic gel. (lid also filled with gel on some models)



055411



Lid	Transparent PC	Transparent PC	PC white	Transparent PC	PC white	Transparent PC	Transparent PC
Body colour	Blue	Blue	Blue	Green	Green	Green	Green
Gel lid	No	No	Yes	No	Yes	No	No
Number of tubes	3 x 4	4 x 5	4 x 8	3 x 4	4 x 8	3 x 4	3 x 4
L x W x H (mm)	151 x 108 x 125	225 x 165 x 105	243 x 157 x 146	151 x 108 x 125	243 x 157 x 146	197 x 140 x 190	197 x 140 x 190
Tube type	Microtubes 0,5/ 1,5 ml	Cryotubes, microtubes	Microtubes 0,5/ 1,5 ml	Microtubes 0,5/ 1,5 ml	Microtubes 0,5/ 1,5 ml	12 / 13 mm tubes	16 / 17 mm tubes
Temperature	-15 °C, 2 hours	-15 °C, 2 hours	-15 °C, 2 hours	below 1°C, 3 hours	below 1°C, 5 hours	below 1°C, 5 hours	below 1°C, 5 hours
Cat. No.	028057	055411	028058	028059	028060	028061	028062
€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	NC -	NC -