

Stericup® Complete filter units

Membrane Express Plus: high flow and low adsorption membrane.

Vacuum filter units for volumes from 150 ml to 1 l.

- Sterile, individual packaging
- Support: polystyrene
- Filtration area : 40 cm²
- Operating temperature : 45 °C maximum
- Single use
- Membranes
 - PES (GP): Express Plus fast flow membrane
 - Durapore (GV): low protein adsorption membrane



Volume receptacle (ml)	Volume funnel (ml)		Stericup VP	Stericup GP	Stericup GV	Stericup HV
		Membrane	Express Plus Membrane (PES) 0.1 µm	Express Plus Membrane (PES) 0.22 µm	Durapore Membrane (PVDF) 0.22 µm	Durapore Membrane (PVDF) 0.45 µm
150	150	Cat. No.	-	051244B	051249B	051284B
		€/box of 12	-	NC -	NC -	NC -
250	250	Cat. No.	051243	051245B	051250B	051285B
		€/box of 12	NC -	NC -	NC -	NC -
500	500	Cat. No.	-	051246B	051251B	051286B
		€/box of 12	-	NC -	NC -	NC -
1000	500	Cat. No.	051955B	051247B	051252B	-
		€/box of 12	NC -	NC -	NC -	-
1100	1000	Cat. No.	-	051248B	051253B	051288B
		€/box of 12	-	NC -	NC -	NC -

Pre-filter for difficult solutions

Cat. No.	Description	Units/box	€/box
051259	Pre-filter for difficult solution	100 membranes	NC -

Stericup E and Steritop E complete filter units

- Membrane Express plus: high flow and low adsorption membrane
- Sterile, individual packaging
- Membrane : 0.22 µm, PES
- Support: polystyrene
- Filtration area : 40 cm²
- Operating temperature : 50 °C
- Single use
- Direct attachment to 38 or 45 mm neck bottles, elimination of funnel reduces waste
- Stericup: filter unit and storage bottle
- Steritop: filter unit, 45 mm neck storage bottle to be ordered separately

Model	Receptacle volume (ml)	Thread (mm)	Cat. No.	€
Stericup	500	38	236055	NC -
		45	236056	NC -
	1000	38	236057	NC -
		45	236058	NC -
Steritop	2,000	38	236059	NC -
		45	236060	NC -

Accessories

CatNo.	Description	€ 12 units
051273B	Empty sterile PS bottle with 45 mm neck without filter 500 ml	NC -
051274B	Empty sterile PS bottle with 45 mm neck without filter 1000 ml	NC -



FILTRATION